

Which Chromatrap® kit is for you?

Are you working with insect samples?

Yes

No

Do you need NGS library reagents?

Are you fixing your samples?

Yes

No

Yes

No

Chromatrap
Drosophila
UniqSeq
(500276/500277)
ProA/G

Chromatrap
Drosophila
ChIP-Seq
(500274/500275)
ProA/G

Are you using FFPE tissue?

Yes

No

Chromatrap
Native
ChIP-Seq
(500237/500238)
ProA/G

Chromatrap FFPE
ChIP-Seq
(500235/500236)
ProA/G

Do you need NGS library reagents?

Yes

No

Are you using Sonication or
Enzymatic digestion?

Sonication

Enzymatic

Are you using Sonication or
Enzymatic digestion?

Sonication

Enzymatic

Chromatrap
UniqSeq
Sonication
(500264/500265)
ProA/G

Chromatrap
UniqSeq
Enzymatic
(500266/500267)
ProA/G

Chromatrap
ChIP-Seq
Sonication
(500189/500190)
ProA/G

Chromatrap
ChIP-Seq
Enzymatic
(500191/500192)
ProA/G



Chromatrap

www.chromatrap.com